



2718 Industrial Drive  
Junction City, KS 66441  
Tel: 800-916-8311  
Fax: 785-238-5688  
[www.InVitria.com](http://www.InVitria.com)  
[info@InVitria.com](mailto:info@InVitria.com)

## Albumin-DX (Lipid-Reduced)

### Product Specifications

**Product Name:** Albumin – DX (Lipid-Reduced)(Recombinant Albumin)  
**Product Number:** 777HSA027  
**Product Grade:** ISO 9001  
**Product Form:** Powder  
**Country of Origin:** U.S.A.  
**Species:** *Oryza sativa* (rice)  
**Tissue:** Seed

No animal components are used in the manufacture of the product, and it is therefore considered as an animal-free material.

---

TEST	SPECIFICATION
Purity (SDS-PAGE)	> 90%
pH (1% solution)	pH 6.5 – 8.5
Solubility (1% solution)	< 0.5 (A600nm)
Conductivity (1% solution)	< 4 mS

---