



2718 Industrial Drive
Junction City, KS 66441
Tel: 800-916-8311
Fax: 785-238-5688
www.InVitria.com
info@InVitria.com



Cellastim – cGMP

Product Specifications

Product Name: Cellastim™ - cGMP (Recombinant Human Serum Albumin)
Product Number: 777HSA024
Product Grade: cGMP
Product Form: Powder
Country of Origin: U.S.A.
Species: *Oryza sativa* (rice)
Tissue: Seed

No animal components are used in the manufacture of the product, and it is therefore considered as an animal-free material.

TEST	SPECIFICATION
PURITY (SEC-HPLC)	≥ 95%
ENDOTOXIN	≤ 4 EU/mg
pH (10% solution)	pH 6.5 – 8.5
SOLUBILITY (10% solution)	< 0.5 (A600nm)
