



2718 Industrial Drive
Junction City, KS 66441
Tel: 800-916-8311
Fax: 785-238-5688
www.InVitria.com
info@InVitria.com

Lacromin™

Product Specifications

Product Name: Lacromin™ (Recombinant Human Lactoferrin - Holo)
Product Number: 777LAC015
Country of Origin: USA
Species: *Oryza sativa* (rice)
Tissue: Seed

No animal components are used in the manufacture of the product, and it is therefore considered as an animal-free material

TEST	SPECIFICATION
Appearance	Pink powder
Lactoferrin content	> 85%
Purity (SDS-PAGE)	≥ 90%
Moisture content	< 10%
pH (1% solution)	6.0 – 8.0
Conductivity (1% solution)	< 4 mS
Bound iron	> 1.2 mg/g LF
Endotoxin	< 10 EU/mg
