



2718 Industrial Drive
Junction City, KS 66441
Tel: 800-916-8311
Fax: 785-238-5688
www.InVitria.com
info@InVitria.com



Product Specifications

Product Name: Cellastim™ (Recombinant Human Serum Albumin)
Product Number: 777HSA017
Product Grade: Cell culture
Product Form: Powder
Country of Origin: USA
Species: *Oryza sativa* (rice)
Tissue: Seed

No animal components are used in the manufacture of the product, and it is therefore considered as an animal-free material.

TEST	SPECIFICATION
PURITY (SDS-PAGE)	> 96%
TOTAL PROTEIN (Nitrogen Analyzer)	≥ 90%
ENDOTOXIN (Pyrogene rFC)	≤ 1 EU/mg
pH (10% solution)	pH 6.5 – 7.5
SOLUBILITY (10% solution)	< 0.5 (A600nm)
