



2718 Industrial Drive
Junction City, KS 66441
Tel: 800-916-8311
Fax: 785-238-5688
www.InVitria.com
info@InVitria.com



Product Specifications

Product Name: Lysobac™ (Recombinant Human Lysozyme)
Product Number: 777LYS016
Product Form: Powder
Country of Origin: USA
Species: *Oryza sativa* (rice)
Tissue: Seed

No animal components are used in the manufacture of the product, and it is therefore considered as an animal-free material

TEST	SPECIFICATION
Protein content ($A_{280nm} 1\% = 24.6$)	$\geq 75\%$
Specific Activity	$> 100,000$ Units/mg protein
Moisture content	$< 10\%$
Solubility (1%, A_{600nm})	< 2.5
