




2718 Industrial Drive
Junction City, KS 66441
Tel: 800-916-8311
Fax: 785-238-5688
www.InVitria.com
info@InVitria.com



Product Specifications

Product Name: Optiferrin™ (Recombinant Transferrin)
Product Number: 777TRF029
Form: Powder
Country of Origin: USA
Species: *Oryza sativa* (rice)
Tissue: Seed

No animal components are used in the manufacture of the product, and it is therefore considered as an animal-free material

TEST	SPECIFICATION
Protein content (Nitrogen analyzer)	> 90%
Purity (SDS-PAGE)	> 95%
Endotoxin	≤ 1 EU/mg
pH (1% solution)	6.0 – 8.0
Solubility (1%, A600 _{nm})	< 0.5
Moisture Content	< 10%
