



2718 Industrial Drive
Junction City, KS 66441
Tel: 800-916-8311
Fax: 785-238-5688
www.InVitria.com
info@InVitria.com

ZAP-CHO™

Product Specifications

Product Name: ZAP-CHO™
Product Number: 777ZAP023
Country of Origin: USA
Species: Non-animal

No animal components are used in the manufacture of the product, and it is therefore considered as an animal-free material.

TEST	SPECIFICATION
Total Protein	≥ 85%
Endotoxin	< 2 EU/mg
pH (1% solution)	pH 6.0 – 8.0
Conductivity (1% solution)	< 4mS
Moisture	< 10%
